


## ELECTROMYOGRAM (EMG)

<https://www.hopkinsmedicine.org/electromyography>

**FUNCTION:** Measures muscle response or electrical activity in response to a nerve's stimulation of the muscle. The test is used to help monitor nerve activity or detect neuromuscular abnormalities.

### TYPES of EMG

- Surface electrodes: They are used for nerve conduction velocity test, to detect true nerve disorders or conditions whereby muscles are affected by nerve injury.



- Needle electrodes: The electrical activity is detected by a needle which serves as an electrode. A needle is inserted through the skin into the muscle.



### **USER INTERFACE**

The activity is displayed visually on a computer screen and may also be detected audibly with a speaker. The presence, size, and shape of the wave form produced on the oscilloscope provide information about the ability of the muscle to respond to nervous stimulation.

### **MAIN ELEMENTS**

- Patients
- Doctors
- Specialists

### **COMPONENTS/DEVICES**

- Computer Screen
- Electrodes
- Microcontrollers
- Amplifiers
- Sensors
- Bluetooth (communication protocol)

### **REFERENCE**

- [https://www.hopkinsmedicine.org/health/treatment-tests-and-therapies/electromyography-emg#vm\\_A\\_cb2f0571](https://www.hopkinsmedicine.org/health/treatment-tests-and-therapies/electromyography-emg#vm_A_cb2f0571)